



Features:

- Only bed alarm capable of setting a “point of alarm” when a patient is in a medically compromising position
- 30-Day, 90-Day, & 1-Year sensors
- Automatically begins patient monitoring within four seconds of patient entering the bed
- Allows staff to compensate for the active or sedate patient by selecting 0, 4, or 8 seconds before an alarm is activated
- Allows staff to select “point of alarm” for individual patient needs
- At a glance, day or night staff can view monitor status modes. Includes IN, OUT, ALARM and STAND-BY
- The patient call cord is inserted into the Informer Plus, eliminating messy and expensive adapters